



Normalization of microarray data

2-channel microarrays

Why do we need to normalize ?

- Microarrays can suffer from various experimental artifacts.
 - Spotting effects (block effect): the spot size (and thus DNA concentration) may be affected by the spotting head - > vary between printing blocks.
 - Position effects: some regions of the slide are more intense, both for the signal and for the background.
 - Dust : dust is highly fluorescent.
 - Dirt: the slide hereby is clearly contains some dirty zone.
 - Scratches: scratches on the glass create autofluorescent zones.
 - Undetectable spots: for some spots, the signal is so small that it can be lower than the surrounding background.
 - Intensity effect: the $\text{Intensity} = f([\text{RNA}])$ curve is not perfectly linear, and is different for the green and red dyes.






Source: SGD, Ogawa, chip : 1562

The raw measurements

Raw measurements

- Raw measurements are provided as mean and background intensities for the red and green channels.
- For each spot on the slide :
 - $R = R_m - R_b$
 - R_m red mean
 - R_b red background
 - $G = G_m - G_b$
 - G_m green mean
 - G_b green background

Never use ratios



- $r=R/G$
- The ratio is a very poor statistics.
- It reflects very badly the regulation :
 - A 10-fold up-regulation is represented by a value of 10, its distance to the random expectation is 9.
 - A 10-fold down-regulation is represented by a ratio of 0.1. Its distance to the random expectation is 0.9.
- Raw ratios will thus emphasize up-regulation, and ignore down-regulation.

Log-ratios

- Using log-ratios has a normalizing effect.
- Usually, a base 2 is used for the log, because it is more intuitive and easy to convert.
- Some examples:


□ R & G	ratio	log(ratio)	regulation
□ R=G	1	0	random expectation
□ R=G*2	2	1	2-fold up-regulation
□ R=G*4	4	2	4-fold up-regulation
□ R=G/4	0.25	-2	4-fold down-regulation
- The statistic is symmetric: up- and down-regulated genes are at the same distance from random expectation (0)

Log ratio distribution

- Channel bias


This chip is visibly not centred around zero. The negative trends suggests a bias towards green channel.

Log ratio distribution; chip 3




Biases in microarray samples

RG (Red-Green) plot




- Channel bias
 - The majority of the dots are below the diagonal.
- Intensity effect
 - The cloud is curved, suggesting a non-linear response of red and/or green channels.
 - A linear regression does not fit well the cloud.

MA plot = RI (Ratio-Intensity) plot




- M is the log-ratio
 - $M = \log(R/G)$
- A is the average log intensity
 - $A = \log(R*G)/2$
 $= [\log(R) + \log(G)]/2$
- The MA plot emphasizes the bad centring and the intensity bias.
- Channel bias
 - The mean ratio and median log-ratio differ from 0.
- Intensity effect
 - The cloud is visibly curved
 - The regression line does not fit the cloud.

Locally weighted linear regression (lowess)



- Locally weighted regression (**LOWESS**) consists in calculating, for each X value, the regression line on the basis of a subset of points around this X value.


Lowess - smooth factor



- The main parameter for lowess is the **smooth factor**, which gives the proportion of points in the plot which influence the smooth at each value. Smaller smooth factor values give a closer fit.
- When the smooth factor is 1, all points influence the regression at each value of the X axis.
- The regression is however not linear, because the influence of a point on another one diminishes with the distance (gaussian kernel).


Lowess - smooth factor

- A smooth factor of 0.5 fits quite well the curve.




Lowess - smooth factor

- With a smooth factor of 0.1, the regression line shows irregularities.



Lowess - smooth factor




- With a smooth factor of 0.01, the regression line is very irregular.
- At each position, the regression is calculated with a small number of neighbours, and is thus strongly influenced by local fluctuations.
- The curve follows local fluctuations of the cloud, which are likely to reflect random effects rather than intrinsic variations.
- This is a problem of **over-fitting**.

Normalization methods


- The BioConductor library `sma` contains various methods of normalization.
 - Median centring
 - Global (chip-wise) LOWESS
 - Block-wise LOWESS
 - Block-wise LOWESS with scaling
- The function `plot.mva()` performs either of these normalizations and draws the MA plot of the normalized data.

Without normalization



- In the next slides, we will apply the different normalization methods to the same chip, and comment the result.
- The original chip is
 - biased towards the green channel
 - Intensity-biased


Median centring



- The median of the log-ratios is subtracted from each log-ratio value.

$$M_{norm} = M - \text{median}(M)$$


Global (chip-wise) lowess normalization



- For chip-wise LOWESS, a regression curve $y(A)$ is calculated with all the spots of the chip.
- The values are normalized by subtracting the regression curve from the M value.

$$M_{norm} = M - y(A)$$


Block-wise lowess (one lowess per print tip)



- For block-wise LOWESS, a regression curve $y_{block}(A)$ is calculated for each block separately.
- The values are normalized by subtracting the block-specific regression curve from the M value.

$$M_{norm} = M - y_{block}(A)$$


Block-wise lowess with scaling



- The block-wise LOWESS can be combined to a scaling operation: each value is divided by a block-specific estimator of dispersion s_{block}
- Yang et al. (Technical report 589) recommend to use the median absolute deviation (MAD) to estimate the block-specific dispersion.

$$M_{norm} = \frac{M - y_{block}(A)}{s_{block}(A)}$$


One lowess curve per print tip




- BloConductor contains a function `plot.print.tip.lowess()` which draws a MA plot with one regression curve per print tip
- This function is convenient to check the result of the normalization, and to compare the different normalization methods.

Comparison between normalization methods


Without normalization




Median centering




Global lowess



Print tip lowess



Scaled print tip lowess



***Chip-wise selection of significant genes
in 2-channel microarrays***

Filtering genes on the basis of their log-ratio


- In the first publications on microarray analysis, genes were filtered on the basis of a threshold on the log-ratio. Typically, papers from Stanford were considering as significantly regulated all genes with
 - R/G $\log_2(R/G)$ regulation
 - ≥ 2 ≥ 1 up-regulated
 - $\leq 1/2$ ≤ -1 down-regulated
- These thresholds were based on an empirical observation (a control chip). They however suffer from several drawbacks
 - They do not rely on any statistical or probabilistic criterion.
 - They do not take into account the bias in centring.
 - This can be circumvented by first centring each chip independently.
 - They do not take into account the chip-specific dispersion. Among a series, some chips may have a wider dispersion than others, due to experimental bias (scanner setting, problems with dye, ...).
 - A scaling is thus required, but after scaling, the values do not directly represent expression ratios anymore.

Significance testing

- We can evaluate the significance of each observation, by calculating its P-value.

$$Pvalue = P(X \geq x)$$

- Under the assumption of normality, the P-value can be obtained from z-scores. Z-scores represent the number of standard deviations from the mean.



$$z = (x - m) / s$$
$$Pvalue = P(Z \geq z)$$

Bonferroni rule

- Multi-testing
 - Assessing the significance of each gene on a chip represents thousands of simultaneous tests. Let N be the number of genes.
 - The risk of error (P-value) associated to each gene will thus be challenged N times.
 - The significance thresholds used for single testing (0.01, 0.001) are thus likely to return many false positive.
- Bonferroni rule
 - Adapt the threshold to the number of simultaneous tests.

$$\alpha \leq \frac{1}{N}$$

E-value

- An alternative but equivalent way to treat the problem of multi-testing is to calculate the expected value for each observation.
- One can then select a threshold on E-value according to the number of false positive considered as acceptable.

$$Evalue = Pvalue * N$$

References

- Normalization

- Yang, Y. H., Dudoit, S., Luu, P., Lin, D. M., Peng, V., Ngai, J. & Speed, T. P. (2002). Normalization for cDNA microarray data: a robust composite method addressing single and multiple slide systematic variation. *Nucleic Acids Res* 30(4), e15.

- Data sources

- Callow, M. J., Dudoit, S., Gong, E. L., Speed, T. P. & Rubin, E. M. (2000). Microarray expression profiling identifies genes with altered expression in HDL-deficient mice. *Genome Res* 10(12), 2022-9.